

For Research Use Only

HNRNPD Monoclonal antibody, PBS Only



Catalog Number: 68236-1-PBS

Basic Information

Catalog Number: 68236-1-PBS	GenBank Accession Number: BC026015	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 3184	CloneNo.: 2G9E6
Source: Mouse	UNIPROT ID: Q14103	
Isotype: IgG2b	Full Name: heterogeneous nuclear ribonucleoprotein D (AU-rich element RNA binding protein 1, 37kDa)	
Immunogen Catalog Number: AG3458	Calculated MW: 355 aa, 38 kDa	
	Observed MW: 38 kDa	

Applications

Tested Applications:
WB, IF, IHC, ELISA

Species Specificity:
Human, mouse, rat, pig, rabbit

Background Information

The heterogeneous nuclear ribonucleoprotein (hnRNP) complexes which provide the substrate for the processing events that pre-mRNAs undergo before becoming functional, translatable mRNAs in the cytoplasm. HNRNPD, also known as AUF1, belongs to the family of AU-binding proteins (AU-BPs) that regulate the cellular half-lives of many mRNAs by directly interacting with an AU-rich element (ARE) located in their 3' untranslated region (PMID:8578590,12704645). AUF1 has four isoforms produced by alternative splicing of a single transcript: p37, p40, p42, and p45 (PMID:9521873).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.


Storage Buffer:
PBS only

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 68236-1-Ig (HNRNP D antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 68236-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using 68236-1-Ig (HNRNP D antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 68236-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using 68236-1-Ig (HNRNP D antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 68236-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using HNRNP D antibody (68236-1-Ig, Clone: 2G9E6) at dilution of 1:400 and CoraLite@488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L), CL594-Phalloidin (red). This data was developed using the same antibody clone with 68236-1-PBS in a different storage buffer formulation.