

For Research Use Only

# APC Anti-Human STRO-1 (STRO-1)

Catalog Number: APC-65184



## Basic Information

Catalog Number:

APC-65184

Size:

100tests , 5 µl/test

Source:

Mouse

Isotype:

IgM, lambda

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

The purified antibody is conjugated with allophycocyanin (APC) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

CloneNo.:

STRO-1

Excitation/Emission maxima wavelengths:

650 nm / 660 nm

## Applications

Tested Applications:

FC

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> human PBMCs were surface stained with 5 uL APC Anti-Human STRO-1 (APC-65184, Clone:STRO-1) or Mouse IgM Isotype Control. Cells were not fixed. CD45 dim, CD34 positive lymphocytes were gated.