

For Research Use Only

CoraLite® Plus 647 Anti-Mouse Integrin beta-7 (FIB504)



Catalog Number:CL647-65213

Basic Information

Catalog Number: CL647-65213	GenBank Accession Number: BC011184	Purification Method: Protein A purification
Size: 100ug, 500 µg/ml	GeneID (NCBI): 16421	CloneNo.: FIB504
Source: Rat	UNIPROT ID: P26011	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Isotype: IgG2a, kappa	Full Name: integrin beta 7	

Applications

Tested Applications:
FC

Species Specificity:
Human, Mouse

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 647 Anti-Mouse Integrin beta-7 (CL647-65213, Clone: FIB504) or 0.5 ug CoraLite® Plus 647 Rat IgG2a Isotype Control (2A3) (CL647-65209, Clone: 2A3), and 0.5 ug CoraLite® Plus 488 Anti-Mouse CD3 (17A2) (CL488-65077, Clone: 17A2). Cells were not fixed.