

For Research Use Only

CoraLite® Plus 750 Anti-Human CD86 (BU63)



Catalog Number: **CL750-65165**

Basic Information

Catalog Number: CL750-65165	GenBank Accession Number: BC040261	Purification Method: Protein G purification
Size: 100tests , 5 µl/test	GeneID (NCBI): 942	CloneNo.: BU63
Source: Mouse	Full Name: CD86 molecule	Excitation/Emission maxima wavelengths: 755 nm / 780 nm
Isotype: IgG1, kappa	Calculated MW: 329 aa, 38 kDa	

Applications

Tested Applications:
FC

Species Specificity:
Human

Background Information

CD86 (also known as B7.2) is a costimulatory molecule belonging to the immunoglobulin superfamily. Primarily expressed on antigen-presenting cells (APCs), including B cells, dendritic cells, and macrophages, CD86 is the ligand for two proteins at the cell surface of T cells, CD28 antigen and cytotoxic T-lymphocyte-associated protein 4. Binding of CD86 with CD28 antigen is a costimulatory signal for activation of the T-cell. Binding of CD86 with cytotoxic T-lymphocyte-associated protein 4 negatively regulates T-cell activation and diminishes the immune response.

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 750 Anti-Human CD86 (CL750-65165, Clone:BU63). Cells were not fixed. Monocytes were gated.



1X10⁶ human PBMCs were surface co-stained with FITC Anti-Human CD14 and 5 ul CoraLite® Plus 750 Anti-Human CD86 (CL750-65165, Clone:BU63). Cells were not fixed. Monocytes were gated.