

For Research Use Only

# CoraLite® Plus 750 Anti-Mouse CD94 (18D3)



Catalog Number: CL750-65240

## Basic Information

<b>Catalog Number:</b> CL750-65240	<b>GenBank Accession Number:</b> BC120854	<b>Purification Method:</b> Affinity purification
<b>Size:</b> 100ug, 0.5 mg/ml	<b>GeneID (NCBI):</b> 16643	<b>CloneNo.:</b> 18D3
<b>Source:</b> Rat	<b>Full Name:</b> killer cell lectin-like receptor, subfamily D, member 1	<b>Excitation/Emission maxima wavelengths:</b> 755 nm / 780 nm
<b>Isotype:</b> IgG2a, kappa		

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
Mouse

## Background Information

### Storage


**Storage:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> mouse splenocytes were surface stained with FITC Anti-Mouse CD49b (HMA2) (FITC-65245, Clone: HMA2) and 0.5 ug CoraLite® Plus 750 Anti-Mouse CD94 (CL750-65240, Clone: 18D3) or CoraLite® Plus 750 Rat IgG2a isotype control. Cells were not fixed.