

For Research Use Only

FITC Anti-Mouse CD115 (c-fms) (AFS98)



Catalog Number: **FITC-65095**

Basic Information

Catalog Number: FITC-65095	GenBank Accession Number: BC036343	Purification Method: Affinity purification
Size: 100ug, 0.5 mg/ml	GeneID (NCBI): 12978	CloneNo.: AFS98
Source: Rat	Full Name: colony stimulating factor 1 receptor	Excitation/Emission maxima wavelengths: 494 nm / 520 nm
Isotype: IgG2a, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ Balb/c mouse peritoneal macrophages were surface stained with 0.5 ug FITC Anti-Mouse CD115 (c-fms) (FITC-65095, Clone: AFS98) (green) or 0.5 ug FITC-Rat IgG2a Isotype Control (black). Cells were not fixed.