

For Research Use Only

FITC Anti-Mouse TCR Beta (H57-597)

Catalog Number: FITC-65106



Basic Information

| | | |
|--------------------------------------|---|---|
| Catalog Number: FITC-65106 | GenBank Accession Number: GeneID (NCBI): 21577 | Purification Method: Affinity purification |
| Size: 100ug, 0.5 mg/ml | Full Name: T-cell receptor beta chain | CloneNo.: H57-597 |
| Source: Armenian Hamster | | Excitation/Emission maxima wavelengths: 494 nm / 520 nm |
| Isotype: IgG | | |

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

T-cell receptor (TCR) consists of heterodimeric glycoproteins (TCR alpha/beta or gamma/delta) that demonstrate homology with immunoglobulins (PMID: 2144901). TCR beta is a member of the immunoglobulin superfamily and is expressed on alpha/beta TCR-bearing T cells and thymocytes. The TCR heterodimers are noncovalently associated with CD3. The TCR/CD3 complex of T-lymphocytes consists of either a TCR alpha/beta or TCR gamma/delta heterodimer coexpressed at the cell surface with the invariant subunits of CD3 labeled gamma, delta, epsilon, zeta, and eta (PMID: 1826255). The TCR recognizes antigens bound to major histocompatibility complex (MHC) molecules. TCR-mediated peptide-MHC recognition is transmitted to the CD3 complex, leading to the intracellular signal transduction (PMID: 11985657).

Storage


Storage:
Store at 2-8°C. Avoid exposure to light.

Storage Buffer:
Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ C57BL/6 mouse splenocytes were surface stained with 0.25 ug FITC-Anti-Mouse TCR beta (FITC-65106, clone H57-597) (red) or 0.25 ug FITC-Armenian Hamster IgG isotype control (blue). Cells were not fixed.