

For Research Use Only.  
Not For Use In Diagnostics.

# Rhodamine (TRITC)-conjugated Goat Anti-Human IgG(H+L)

Catalog Number:SA00007-10



## Information

### Catalog Number:

SA00007-10

### Host:

Goat

### Applications:

IF, FC

### Reactivity:

Human

### Physical State:

Liquid

### Conjugation:

Rhodamine

### Clonality:

Polyclonal

## Recommended Dilutions

1:20-1:100 for most applications

## Preservative

0.02% Sodium Azide

## Fluorophore

Tetramethyl Rhodamine Isothiocyanate (TRITC) Amax=550nm, Emax=570nm

## Stabilizer

5 mg/ml BSA

## Absorbance Ratio

A550/A280=0.48

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

