

For Research Use Only.  
Not For Use In Diagnostics.

# CoraLite488-conjugated Mouse Anti-Heavy Chain of Rabbit IgG

Catalog Number:SA00014-4



[www.ptglab.com](http://www.ptglab.com)

## Information

Catalog Number:

SA00014-4

Host:

Mouse

Applications:

IF, FC

Reactivity:

Rabbit

Physical State:

Liquid

Conjugation:

CoraLite®488

## Recommended Dilutions

1:100-1:500

## Preservative

Proclin300

## Fluorophore

CL488 Fluorescent Dye, Amax=488 nm; Emax=515 nm

## Stabilizer

0.5% BSA

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

## Purity

Protein G purification

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data