

For Research Use Only.
Not For Use In Diagnostics.

CoraLite647-conjugated AffiniPure F(ab')₂ Fragment Donkey Anti- Mouse IgG (H+L)



Catalog Number:SA00014-8

Information

Catalog Number: SA00014-8	Reactivity: Mouse
Host: Donkey	Physical State: Liquid
Applications: IF, FC	Conjugation: CoraLite®647

Recommended Dilutions

1:100 – 1:500

Preservative

0.01% Sodium Azide

Fluorophore

CL647 Fluorescent Dye. A_{max}=642nm, E_{max}=662nm

Stabilizer

5 mg/ml BSA

Absorbance Ratio

A₆₅₁/A₂₈₀ = 4.0

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.01% Sodium Azide, pH 7.6.
Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)


E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using 66009-1-Ig (Beta Actin antibody) and CoraLite®647 conjugated AffiniPure F(ab')₂ Fragment Donkey Anti-Mouse IgG (H+L) antibody (SA00014-8) at dilution of 1:150.



1X10⁶ HepG2 cells were intracellularly stained with 0.2 ug Anti-human RBAP48 Monoclonal antibody (66060-1-Ig) and CoraLite®647-Conjugated AffiniPure F(ab')₂ Fragment Donkey Anti-Mouse IgG (H+L) (SA00014-8) at dilution 1:1000 (red), and 0.2 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).