

Allgemeine Informationen

Katalog-Nr.:

10152-2-AP

Größe:

150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG0204

GenBank-Zugangsnummer:

BC002781

GeneID (NCBI):

6429

Vollständiger Name:

splicing factor, arginine/serine-rich 4

Berechnete Masse:

57 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

In Publikationen genannte Anwendungen:

IP, WB

Getestete Reaktivität:

Human, Maus, Ratte

Zitierte Arten:

Human, Maus

Hintergrundinformationen

Splicing factor, arginine/serine-rich 4 (SFRS4, synonym: SRP75). SR proteins are a family of proteins that have a common epitope that binds active sites of polymerase II transcription. SRp75 protein, like other SR proteins, contains an N-terminal RNA recognition motif (RRM), a glycine-rich region, an internal region homologous to the RRM, and a long (315-amino-acid) C-terminal domain composed predominantly of alternating serine and arginine residues. Four of the other SR family members have been shown to restore activity to an otherwise splicing-deficient extract. SRp75 together with SRp20 can also complement the splicing-deficient extract. The apparent molecular mass of dephosphorylated SRp75 is 57 kDa.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
C M Exline	18272582	J Virol	IP
Yuanbiao Guo	19865483	PLoS One	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

