

Nur für Forschungszwecke

TBRG4 Polyklonaler Antikörper

Katalog-Nr.: 10177-2-AP



Allgemeine Informationen

Katalog-Nr.:

10177-2-AP

Größe:

150ul, Konzentration: 2000 µg/ml von 9238

Nanodrop und 620 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG0212

GenBank-Zugangsnummer:

BC002732

GeneID (NCBI):

9238

Vollständiger Name:

transforming growth factor beta regulator 4

Berechnete Masse:

71 kDa

Beobachtete Masse:

58 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : A431-Zellen, HepG2-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



A431 cells were subjected to SDS PAGE followed by western blot with 10177-2-AP (TBRG4 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.