

FMNL1 Polyklonaler Antikörper

Katalog-Nr.:10395-1-AP

Allgemeine Informationen

Katalog-Nr.: 10395-1-AP	GenBank-Zugangsnummer: BC001710	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 245 µg/ml von Nanodrop;	GenID (NCBI): 752	Empfohlene Verdünnungen: WB 1:200-1:500
Wirt: Kaninchen	Vollständiger Name: formin-like 1	
Isotyp: IgG	Berechnete Masse: 160 kDa	
Immunogen Katalognummer: AG0594	Beobachtete Masse: 150-160 kDa	

Anwendungen

Geprüfte Anwendungen: WB,ELISA	Positivkontrollen: WB : Raji-Zellen, HeLa-Zellen, Jurkat-Zellen
Getestete Reaktivität: Human	

Hintergrundinformationen

FMNL1, also named as C17orf1, C17orf1B and FMNL, belongs to the formin homology family. It may play a role in the control of cell motility and survival of macrophages. FMNL protein may be involved in anti-apoptosis mechanisms. It is overexpressed in lymphoid malignancies and is found associated with Akt. The MW of FMNL1 is about 150-160 kDa. (PMID: 16494944)

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Raji cells were subjected to SDS PAGE followed by western blot with 10395-1-AP (FMNL1 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.