

Allgemeine Informationen

Katalog-Nr.:

10451-1-AP

Größe:

150ul , Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG0660

GenBank-Zugangsnummer:

BC002806

GeneID (NCBI):

8612

Vollständiger Name:

phosphatidic acid phosphatase type
2C

Berechnete Masse:

32 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB : Mauslebergewebe,

Hintergrundinformationen

PPAP2C(phosphatidic acid phosphatase type 2C) is also named as PAP2C, LPP2, PAP-2c, PAP2-g. PPAP2C is a member of the phosphatidic acid phosphatase family of enzymes that regulate the dephosphorylation of lipid phosphates and the gene appears to be particularly highly expressed in breast carcinoma, dedifferentiated chondrosarcoma, fibromatosis, synovial sarcoma, and neurofibroma(PMID:19139135). It has 2 isoforms produced by alternative splicing and the full length protein has a glycosylation site.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse liver tissue were subjected to SDS PAGE followed by western blot with 10451-1-AP (PPAP2C Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.