

Nur für Forschungszwecke

HEM1; NCKAP1L Polyklonaler Antikörper

Katalog-Nr.: 10454-1-AP



Allgemeine Informationen

Katalog-Nr.:

10454-1-AP

Größe:

150ul, Konzentration: 247 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG0706

GenBank-Zugangsnummer:

BC001604

GeneID (NCBI):

3071

Vollständiger Name:

NCK-associated protein 1-like

Berechnete Masse:

128 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

Membrane-associated protein HEM1 encoded by human NCKAP1L (Nck-associated protein 1-like) gene belongs to the HEM family of tissue-specific transmembrane proteins which are highly conserved from invertebrates through mammals. HEM1 is characterized as hematopoietic cell-specific actin regulatory protein. HEM1 is a component of the "WAVE" complex, which signals downstream of activated Rac to stimulate F-actin polymerization in response to immuno-receptor signaling. Therefore, HEM1 may be essential for most active processes dependent on reorganization of the actin cytoskeleton in immune cells.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

