

Nur für Forschungszwecke

PRRG2 Polyklonaler Antikörper

Katalog-Nr.: 11618-1-AP



Allgemeine Informationen

Katalog-Nr.:

11618-1-AP

Größe:

150ul, Konzentration: 450 µg/ml von Nanodrop und 373 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG2197

GenBank-Zugangsnummer:

BC026032

GeneID (NCBI):

5639

Vollständiger Name:

proline rich Gla (G-carboxyglutamic acid) 2

Berechnete Masse:

179 aa, 20 kDa

Beobachtete Masse:

16 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB: A549-Zellen,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



A549 cells were subjected to SDS PAGE followed by western blot with 11618-1-AP (PRRG2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.