

Nur für Forschungszwecke

TULP2 Polyklonaler Antikörper

Katalog-Nr.:12668-1-AP



Allgemeine Informationen

Katalog-Nr.:
12668-1-AP

Größe:
150ul, Konzentration: 160 µg/ml von
Nanodrop und 160 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG3350

GenBank-Zugangsnummer:
BC026070

GeneID (NCBI):
7288

Vollständiger Name:
tubby like protein 2

Berechnete Masse:
520 aa, 59 kDa

Beobachtete Masse:
55-59 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung
Empfohlene Verdünnungen:
WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human

Positivkontrollen:

WB : PC-3-Zellen, humanes Lebergewebe, humanes
Lungengewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



PC-3 cells were subjected to SDS PAGE followed by western blot with 12668-1-AP (TULP2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.