

Nur für Forschungszwecke

PAMCI Polyklonaler Antikörper

Katalog-Nr.: 12861-1-AP



Allgemeine Informationen

Katalog-Nr.:

12861-1-AP

Größe:

150ul, Konzentration: 600 µg/ml von
Nanodrop;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG3976

GenBank-Zugangsnummer:

BC031589

GeneID (NCBI):

9182

Vollständiger Name:

Ras association (RalGDS/AF-6)
domain family (N-terminal) member
9

Berechnete Masse:

435 aa, 50 kDa

Beobachtete Masse:

50-55 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:3000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB : HepG2-Zellen, humanes Lebergewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



HepG2 cells were subjected to SDS PAGE followed by western blot with 12861-1-AP (PAMCI antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.