

Nur für Forschungszwecke

KLF7 Polyklonaler Antikörper

Katalog-Nr.: 13125-1-AP

3 Publikationen



Allgemeine Informationen

Katalog-Nr.:

13125-1-AP

Größe:

150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG3797

GenBank-Zugangsnummer:

BC012919

GeneID (NCBI):

8609

Vollständiger Name:

Kruppel-like factor 7 (ubiquitous)

Berechnete Masse:

230 aa, 25 kDa

Beobachtete Masse:

33 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

In Publikationen genannte Anwendungen:

IF, WB

Getestete Reaktivität:

Human, Maus, Ratte

Zitierte Arten:

Human

Hintergrundinformationen

KLF7, also named as UKLF, is a 302 amino acid protein, which contains three C2H2-type zinc fingers and belongs to the krueppel C2H2-type zinc-finger protein family. KLF7 is ubiquitous and highly expressed in brain and spinal cord in the adult, and in kidney and brain in the embryo. KLF7 as a transcriptional activator localizes in the nucleus and binds in vitro to the CACCC motif of the beta-globin promoter and to the SP1 recognition sequence.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Endong Zhu	29348304	J Mol Endocrinol	
Gao Yan Y	22614632	Int J Mol Med	IF
Ning Wang	34012590	J Thorac Dis	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

