

UIMC1 Polyklonaler Antikörper

Katalog-Nr.: 13642-1-AP

Allgemeine Informationen

Katalog-Nr.:

13642-1-AP

Größe:

150ul, Konzentration: 220 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG4582

GenBank-Zugangsnummer:

BC032561

GeneID (NCBI):

51720

Vollständiger Name:

ubiquitin interaction motif containing
1

Berechnete Masse:

719 aa, 80 kDa

Beobachtete Masse:

97 kDa, 80 kDa, 71-75 kDa, 61 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:200-1:1000
IP 0.5-4.0 µg für IP und 1:500-1:1000
für WB
IHC 1:20-1:200
IF 1:10-1:100

Anwendungen

Geprüfte Anwendungen:

IF, IHC, IP, WB, ELISA

Getestete Reaktivität:

Human, Maus

Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.

Positivkontrollen:

WB : Maus-Kolontgewebe, COLO 320-Zellen, HeLa-Zellen, Raji-Zellen

IP : COLO 320-Zellen,

IHC : humanes Hodengewebe,

IF : HeLa-Zellen, MCF-7-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




mouse colon tissue were subjected to SDS PAGE followed by western blot with 13642-1-AP (UIMC1 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.




Immunofluorescent analysis of HeLa cells, using UIMC1 antibody 13642-1-AP at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



IP Result of anti-UIMC1 (IP:13642-1-AP, 4ug; Detection:13642-1-AP 1:500) with COLO 320 cells lysate 2400ug.



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 13642-1-AP (UIMC1 Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 13642-1-AP (UIMC1 Antibody) at dilution of 1:50 (under 40x lens).