

MIB2 Polyklonaler Antikörper

Katalog-Nr.: 13696-1-AP

Allgemeine Informationen

Katalog-Nr.:

13696-1-AP

Größe:

150ul, Konzentration: 227 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG4592

GenBank-Zugangsnummer:

BC037542

GeneID (NCBI):

142678

Vollständiger Name:

mindbomb homolog 2 (Drosophila)

Berechnete Masse:

110 kDa

Beobachtete Masse:

70 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

MIB2, also named as E3 ubiquitin-protein ligase MIB2, is a 1013 amino acid protein, which contains 9 ANK repeats and 2 MIB/HERC2 domains. MIB2 as E3 ubiquitin-protein ligase that mediates ubiquitination of Delta receptors, which act as ligands of Notch proteins. MIB2 positively regulates the Delta-mediated Notch signaling by ubiquitinating the intracellular domain of Delta, leading to endocytosis of Delta receptors. MIB2 is expressed in skeletal muscle, and to a lesser extent in heart, brain and kidney.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

