

Nur für Forschungszwecke

ARHGAP20 Polyklonaler Antikörper

Katalog-Nr.: 13775-1-AP



Allgemeine Informationen

Katalog-Nr.:

13775-1-AP

Größe:

150ul, Konzentration: 187 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG4764

GenBank-Zugangsnummer:

BC039340

GeneID (NCBI):

57569

Vollständiger Name:

Rho GTPase activating protein 20

Berechnete Masse:

133 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : humanes Hirngewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



human brain tissue were subjected to SDS PAGE followed by western blot with 13775-1-AP (ARHGAP20 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.