

Allgemeine Informationen

Katalog-Nr.:

14633-1-AP

Größe:

150ul, Konzentration: 180 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG6212

GenBank-Zugangsnummer:

BC065236

GeneID (NCBI):

23521

Vollständiger Name:

ribosomal protein L13a

Berechnete Masse:

23 kDa

Beobachtete Masse:

25-28 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus

Hintergrundinformationen

Ribosomal protein L13a(RPL13a) is a member of the L13 ribosomal protein family and a structural component of the large 60S ribosomal subunit. It's a reference gene(RG) that stably expressed not only during differentiation, but also during expansion to higher passages. And it's suggested as a RG for a variety research field. RPL13a is a key component of the IFN-gamma-activated inhibitor of translation(GAIT) complex, which inhibits translation of a select group of inflammation-related transcripts via binding to the 3' UTR of target mRNA. once phosphorylated of Ser77 by ZIKV, RPL13a released from the complex.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

