

Nur für Forschungszwecke

ANAPC5 Polyklonaler Antikörper

Katalog-Nr.:14926-1-AP



Allgemeine Informationen

Katalog-Nr.: 14926-1-AP	GenBank-Zugangsnummer: BC001081	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 187 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 51433	Empfohlene Verdünnungen: WB 1:200-1:1000 IHC 1:50-1:500
Wirt: Kaninchen	Vollständiger Name: anaphase promoting complex subunit 5	
Isotyp: IgG	Berechnete Masse: 85 kDa	
Immunogen Katalognummer: AG6688	Beobachtete Masse: 85 kDa	

Anwendungen

Geprüfte Anwendungen: IHC, WB,ELISA	Positivkontrollen:
Getestete Reaktivität: Human, Maus, Ratte	WB : Mauslungengewebe, BxPC-3-Zellen
Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.	IHC : humanes Pankreaskarzinomgewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten



mouse lung tissue were subjected to SDS PAGE followed by western blot with 14926-1-AP (ANAPC5 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 14926-1-AP (ANAPC5 antibody) at dilution of 1:200 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 14926-1-AP (ANAPC5 antibody) at dilution of 1:200 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).