

Nur für Forschungszwecke

PCOLCE Polyklonaler Antikörper

Katalog-Nr.: 14993-1-AP

1 Publikationen



Allgemeine Informationen

Katalog-Nr.:

14993-1-AP

Größe:

150ul, Konzentration: 160 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG6979

GenBank-Zugangsnummer:

BC000574

GeneID (NCBI):

5118

Vollständiger Name:

procollagen C-endopeptidase
enhancer

Berechnete Masse:

48 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

PCOLCE (Procollagen C-endopeptidase enhancer 1), also named as PCPE1, is a secreted glycoprotein that functions as a positive regulator of procollagen processing. PCOLCE binds to the C-terminal propeptide of type I procollagen and enhances procollagen C-proteinase activity. Procollagen C-proteinases have important roles in developmental processes and assembly of the extracellular matrix.

Bemerkenswerte Veröffentlichungen

Verfasser

Pubmed ID

Journal

Anwendung

Efrat Kessler

30840244

Methods Mol Biol

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

