

Nur für Forschungszwecke

SEC61alpha Polyklonaler Antikörper

Katalog-Nr.:15512-1-AP



Allgemeine Informationen

Katalog-Nr.: 15512-1-AP	GenBank-Zugangsnummer: BC002951	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 140 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 29927	
Wirt: Kaninchen	Vollständiger Name: Sec61 alpha 1 subunit (S. cerevisiae)	
Isotyp: IgG	Berechnete Masse: 52 kDa	
Immunogen Katalognummer: AG7693		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

The SEC61 complex is a major component of the endoplasmic reticulum (ER) translocation site (PMID: 7720858). This complex is composed of three subunits: SEC61alpha, SEC61beta and SEC61gamma. SEC61alpha is a polytopic integral membrane protein. It is a core component of the ER translocation complex, playing a crucial role in the insertion of secretory and membrane polypeptides into the ER (PMID: 9512475).

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

