

Nur für Forschungszwecke

SNX15 Polyklonaler Antikörper

Katalog-Nr.:16049-1-AP



Allgemeine Informationen

Katalog-Nr.:
16049-1-AP

Größe:
150ul, Konzentration: 300 µg/ml von
Nanodrop und 213 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG9006

GenBank-Zugangsnummer:
BC009897

GeneID (NCBI):
29907

Vollständiger Name:
sorting nexin 15

Berechnete Masse:
342 aa, 38 kDa

Beobachtete Masse:
50 kDa and 38 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:3000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB : HeLa-Zellen, COLO 320-Zellen, humanes
Hodengewebe, Maushirngewebe, Mauslungengewebe,
SH-SY5Y-Zellen

Hintergrundinformationen

SNX15 belongs to the sorting nexin family and is involved in several stages of intracellular trafficking. Overexpression of SNX15 results in a decrease in the processing of insulin and hepatocyte growth factor receptors to their mature subunits (PMID: 11085978). SNX15 is widely expressed, with the highest expression in skeletal muscle, heart, brain, kidney, spleen, thymus, and small intestine.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 16049-1-AP (SNX15 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.