

Nur für Forschungszwecke

MDM1-Specific Polyklonaler Antikörper

Katalog-Nr.: 16171-1-AP



Allgemeine Informationen

Katalog-Nr.:
16171-1-AP

Größe:

150ul, Konzentration: 700 µg/ml von Nanodrop und 360 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_020128

GeneID (NCBI):

56890

Vollständiger Name:

Mdm1 nuclear protein homolog (mouse)

Berechnete Masse:

81 kDa

Beobachtete Masse:

81 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB : humanes Hirngewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



human brain tissue were subjected to SDS PAGE followed by western blot with 16171-1-AP (MDM1-Specific antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.