

Nur für Forschungszwecke

MAGEC1-Specific Polyklonaler Antikörper

Katalog-Nr.: 16176-1-AP



Allgemeine Informationen

Katalog-Nr.:

16176-1-AP

Größe:

150ul, Konzentration: 473 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_005462

GeneID (NCBI):

9947

Vollständiger Name:

melanoma antigen family C, 1

Berechnete Masse:

124 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

This gene is a member of the melanoma antigen gene (MAGE) family. The proteins of this family are tumor-specific antigens. Expression of this protein was not seen in normal tissues, with the exception of testis. This protein contains a large number of unique short repetitive sequences in front of the MAGE-homologous sequence, and therefore is about 800 aa longer than the other MAGE proteins. (From OMIM)

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

