

Nur für Forschungszwecke

ADAM15-Specific Polyklonaler Antikörper

Katalog-Nr.: 16548-1-AP



Allgemeine Informationen

Katalog-Nr.:

16548-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_207197

GeneID (NCBI):

8751

Vollständiger Name:

ADAM metalloproteinase domain 15

Berechnete Masse:

93 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

ADAM15, a member of the ADAM (a disintegrin and metalloprotease) family, is a membrane protein containing both protease and adhesion domains and may, thus, be involved in tumor invasion and metastasis (PMID:15756594). It is also named as MDC15. ADAM15 is processed during epididymal maturation and acrosome reaction and that it may play a role during sperm-egg binding. Jurkat and testicular cells express the precursor ADAM15 of 110 kDa, as well as the ADAM15 processing form, consisting of about 75-85 kDa (PMID:18390692). It has 11 isoforms produced by alternative splicing. This antibody is specific to ADAM15.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

