

KDEL C1 Polyklonaler Antikörper

Katalog-Nr.: 17018-1-AP

Allgemeine Informationen

Katalog-Nr.:

17018-1-AP

Größe:

150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG9999

GenBank-Zugangsnummer:

BC001297

GeneID (NCBI):

79070

Vollständiger Name:

KDEL (Lys-Asp-Glu-Leu) containing 1

Berechnete Masse:

58 kDa

Beobachtete Masse:

50 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

KDEL C1, also named as EP58 and ERp58, belongs to the KDEL C family. Soluble proteins in the endoplasmic reticulum (ER) contain a specific carboxyterminal sequence KDEL (Lys-Asp-Glu-Leu) and include the coat proteins required for vesicle budding from the ER. KDEL C1 shares no significant homology to any of the known ER-resident proteins. It is primarily expressed in embryo, placenta and adult heart.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with 17018-1-AP (KDEL1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.