

MMP16 Polyklonaler Antikörper

Katalog-Nr.: 17839-1-AP

Allgemeine Informationen

Katalog-Nr.:

17839-1-AP

Größe:

150ul, Konzentration: 193 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG12186

GenBank-Zugangsnummer:

BC075005

GeneID (NCBI):

4325

Vollständiger Name:

matrix metalloproteinase 16
(membrane-inserted)

Berechnete Masse:

607 aa, 70 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte, Zebrafisch

Hintergrundinformationen

MMP16(Matrix metalloproteinase-16) is also named as MMPX2,MT3-MMP.It is an endopeptidase that degrades various components of the extracellular matrix, such as collagen type III and fibronectin.Northern blotting shows MMP16 is expressed as a 12-kb transcript in brain, placenta, heart, and some carcinoma cell lines, but is not detectably expressed in lung, kidney, liver, spleen, and muscle(PMID:7559440).And MMP16 is encountered in endothelial cells of proliferative endometrium in vivo(PMID:15531549).

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

