

Nur für Forschungszwecke

PRKY Polyklonaler Antikörper

Katalog-Nr.:17985-1-AP



Allgemeine Informationen

Katalog-Nr.: 17985-1-AP	GenBank-Zugangsnummer: BC074852	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 207 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 5616	
Wirt: Kaninchen	Vollständiger Name: protein kinase, Y-linked	
Isotyp: IgG	Berechnete Masse: 277 aa, 32 kDa	
Immunogen Katalognummer: AG12404	Beobachtete Masse: 30 kDa	

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

PRKY is a member of the cAMP-dependent serine threonine protein kinase gene family, which appears to encode a protein related to the catalytic subunit of the cAMP-dependent protein kinases, which are key players in the cellular responses to the second messenger cAMP. A chromosomal aberration involving PRKY is a cause of sex reversal disorder. As a homolog with PRKY, the antibody can recognize endogenous level of PRKY and PRKY.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

