

Nur für Forschungszwecke

CYP2R1 Polyklonaler Antikörper

Katalog-Nr.:18180-1-AP



Allgemeine Informationen

Katalog-Nr.:
18180-1-AP

Größe:
150ul, Konzentration: 353 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG12873

GenBank-Zugangsnummer:
BC104909

GeneID (NCBI):
120227

Vollständiger Name:
cytochrome P450, family 2, subfamily
R, polypeptide 1

Berechnete Masse:
501 aa, 57 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



**WB result of CYP2Ra antibody (18180-1-AP, 1:500)
with transfected human CYP2R1 and mouse CYP2R1
by Dr. Jeff Roizen.**