

Nur für Forschungszwecke

ZBTB7A Polyklonaler Antikörper

Katalog-Nr.: 18700-1-AP

1 Publikationen



Allgemeine Informationen

Katalog-Nr.:

18700-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

BC113511

GeneID (NCBI):

51341

Vollständiger Name:

zinc finger and BTB domain
containing 7A

Berechnete Masse:

61 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human, Maus, Ratte

Zitierte Arten:

Human, Maus

Hintergrundinformationen

ZBTB7A, also named as LRF, FBI1, FBI-1, ZBTB7, ZNF857A, pokemon and DKFZp547O146, plays a key role in the instruction of early lymphoid progenitors to develop into B lineage by repressing T-cell instructive Notch signals. ZBTB7A may be involved in transcriptional regulation. It enhances transcription of the NF-kappaB-responsive E-selectin gene by nuclear localization of the p65 subunit of NF-kappaB. ZBTB7A is widely expressed in normal thymus and it is aberrantly overexpressed in human cancers. Results suggest that this protein can specifically stimulate Tat-activated HIV-1 transcription.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Jing-Ping Zhou	31998789	Biomed Res Int	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

