

Nur für Forschungszwecke

LAMC2 Polyklonaler Antikörper

Katalog-Nr.: 19698-1-AP

3 Publikationen



Allgemeine Informationen

Katalog-Nr.:

19698-1-AP

Größe:

150ul, Konzentration: 900 µg/ml von Nanodrop und 413 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_005562

GeneID (NCBI):

3918

Vollständiger Name:

laminin, gamma 2

Berechnete Masse:

131 kDa

Beobachtete Masse:

140 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

In Publikationen genannte Anwendungen:

IHC, WB

Getestete Reaktivität:

Human, Maus, Ratte

Zitierte Arten:

Human, Ratte

Positivkontrollen:

WB: A431-Zellen,

Hintergrundinformationen

LAMC2, also named as LAMB2T and LAMNB2, is thought to mediate the attachment, migration and organization of cells into tissues during embryonic development by interacting with other extracellular matrix components. LAMC2 binds to cells via a high affinity receptor. Laminin exerts cell-scattering activity toward a wide variety of cells, including epithelial, endothelial, and fibroblastic cells.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Giricea Lorusso	36070364	Sci Transl Med	WB, IHC
Caiqun Huang	33179081	Mol Med Rep	WB
Caiming Xu	32244009	J Proteomics	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



A431 cells were subjected to SDS PAGE followed by western blot with 19698-1-AP (LAMC2 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.