

Allgemeine Informationen

Katalog-Nr.:

19893-1-AP

Größe:

150ul , Konzentration: 260 µg/ml von Nanodrop und 200 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG13728

GenBank-Zugangsnummer:

BC002574

GeneID (NCBI):

55159

Vollständiger Name:

ring finger and WD repeat domain 3

Berechnete Masse:

774 aa, 85 kDa

Beobachtete Masse:

85 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

IHC 1:100-1:200

IF 1:20-1:100

Anwendungen

Geprüfte Anwendungen:

IF, IHC, WB, ELISA

In Publikationen genannte Anwendungen:

IF, IHC, WB

Getestete Reaktivität:

Human, Maus, Ratte

Zitierte Arten:

Human

Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.

Positivkontrollen:

WB : Mauslungengewebe,

IHC : humanes Hodengewebe, humanes Kolongewebe

IF : HEK-293-Zellen, HeLa-Zellen

Hintergrundinformationen

RFWD3, also named as RNF201, is a E3 ubiquitin-protein ligase that mediates the ubiquitination of p53/TP53 in the late response to DNA damage, and acts as a positive regulator of p53/TP53 stability, thereby regulating the G1/S DNA damage checkpoint. RFWD3 may act by catalyzing the formation of short polyubiquitin chains on p53/TP53 that are not targeted to the proteasome. In response to ionizing radiation, It interacts with MDM2 and enhances p53/TP53 ubiquitination, possibly by restricting MDM2 from extending polyubiquitin chains on ubiquitinated p53/TP53.(PMID: 20173098) From BLAST result, this antibody is specific to RFWD3. 45kd and 55kd maybe spliceosome of RFWD3.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Yanfei Zhang	31571161	Front Med	IHC
Guochao Liu	33842359	Front Oncol	WB,IF
Ming Gao	33704464	Nucleic Acids Res	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
 W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




mouse lung tissue were subjected to SDS PAGE followed by western blot with 19893-1-AP (RFWD3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human testis using 19893-1-AP (RFWD3 antibody) at dilution of 1:400 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis using 19893-1-AP (RFWD3 antibody) at dilution of 1:400 (under 40x lens).



Immunofluorescent analysis of 293 cells, using RFWD3 antibody 19893-1-AP at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).