

Nur für Forschungszwecke

NT5C1B-Specific Polyklonaler Antikörper

Katalog-Nr.: 20076-1-AP



Allgemeine Informationen

Katalog-Nr.:
20076-1-AP

Größe:
150ul, Konzentration: 287 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

GenBank-Zugangsnummer:
NM_001002006

GeneID (NCBI):
93034

Vollständiger Name:
5'-nucleotidase, cytosolic IB

Berechnete Masse:
69 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

NT5C1B, also named as AIRP, FKSG85, CN-IB, CN1B and MGC26640, belongs to the 5'-nucleotidase type 3 family. NT5C1B dephosphorylates the 5' and 2'(3')-phosphates of deoxyribonucleotides. It helps to regulate adenosine levels. NT5C1B catalyzes the reaction: A 5'-ribonucleotide + H₂O = a ribonucleoside + phosphate. This antibody is specific to NT5C1B. It recognizes all the isoforms of NT5C1B.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

