

Nur für Forschungszwecke

DOCK180-Specific Polyklonaler Antikörper



Katalog-Nr.:20244-1-AP 1 Publikationen

Allgemeine Informationen

Katalog-Nr.: 20244-1-AP	GenBank-Zugangsnummer: NM_001380	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul	GeneID (NCBI): 1793	
Wirt: Kaninchen	Vollständiger Name: dedicator of cytokinesis 1	
Isotyp: IgG	Berechnete Masse: 215 kDa	

Anwendungen

Geprüfte Anwendungen:
ELISA
In Publikationen genannte Anwendungen:
IF
Getestete Reaktivität:
Human
Zitierte Arten:
Human

Hintergrundinformationen

DOCK1, also named as DOCK180, belongs to the DOCK family. It is involved in cytoskeletal rearrangements required for phagocytosis of apoptotic cells and cell motility. DOCK1 functions as a guanine nucleotide exchange factor (GEF), which activates Rac/Rho small GTPases by exchanging bound GDP for free GTP. Its GEF activity may be enhanced by ELMO1. The antibody is specific to DOCK1.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Joseph H R Hetmanski	31607653	Dev Cell	IF

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

