

Nur für Forschungszwecke

# OATP1B3 Polyklonaler Antikörper

Katalog-Nr.:22815-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
22815-1-AP

Größe:  
150ul , Konzentration: 320 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

Immunogen Katalognummer:  
AG18940

GenBank-Zugangsnummer:  
BC105597

GeneID (NCBI):  
28234

Vollständiger Name:  
solute carrier organic anion  
transporter family, member 1B3

Berechnete Masse:  
702 aa, 77 kDa

Reinigungsmethode:  
Antigen-affinitätsgereinigt

## Anwendungen

Geprüfte Anwendungen:  
ELISA

Getestete Reaktivität:  
Human

## Hintergrundinformationen

SLCO1B3, also known as OATP8 or OATP1B3, is a member of Organic Anion Transporting Polypeptides (OATPs) which are sodium-independent organic anion transporters mediating the uptake of a wide range of structurally diverse endogenous and exogenous compounds including bile acids, hormone conjugates, peptides, toxins, as well as a multitude of therapeutic drugs. SLCO1B3 is specifically expressed in liver under normal condition, while its expression has also been observed in cancer tissues like colon, prostate and pancreas.

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

