

ADAM33 Polyklonaler Antikörper

Katalog-Nr.: 22860-1-AP

Allgemeine Informationen

Katalog-Nr.:

22860-1-AP

Größe:

150ul, Konzentration: 187 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG18852

GenBank-Zugangsnummer:

BC125113

GeneID (NCBI):

80332

Vollständiger Name:

ADAM metalloproteinase domain 33

Berechnete Masse:

813 aa, 88 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB: A549-Zellen,

Hintergrundinformationen

ADAM33 is a member of the 'a disintegrin and metalloprotease domain' family of proteins. Endogenous ADAM33 is detected in bronchus tissue, bronchial smooth muscle cells, and MRC-5 fibroblasts, consistent with a role in the pathophysiology of asthma (PMID:12535637). ADAM33 is a novel tumour suppressor gene that may be useful as a molecular marker for invasive lobular carcinoma of the breast (PMID:19267929)

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

