

Nur für Forschungszwecke

ZSCAN21 Polyklonaler Antikörper

Katalog-Nr.: 23164-1-AP



Allgemeine Informationen

Katalog-Nr.:
23164-1-AP

Größe:

150ul, Konzentration: 2000 µg/ml von 7589

Nanodrop und 813 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG19595

GenBank-Zugangsnummer:

BC047309

GeneID (NCBI):

Vollständiger Name:
zinc finger and SCAN domain
containing 21

Berechnete Masse:

473 aa, 54 kDa

Beobachtete Masse:

54 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : HEK-293-Zellen, HepG2-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



HEK-293 cells were subjected to SDS PAGE followed by western blot with 23164-1-AP (ZSCAN21 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.