

BAHD1 Polyklonaler Antikörper

Katalog-Nr.: 23991-1-AP

Allgemeine Informationen

Katalog-Nr.:
23991-1-AP

Größe:
150ul, Konzentration: 400 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG21164

GenBank-Zugangsnummer:
BC022782

GeneID (NCBI):
22893

Vollständiger Name:
bromo adjacent homology domain
containing 1

Berechnete Masse:
780 aa, 85 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

Gene repression in eukaryotes in response to developmental or environmental signals is a complex multistep process requiring the concerted action of many cellular factors, including DNA-binding transcription factors and cofactors. It also involves chromatin components and chromatin-binding/-modifying proteins, which are implicated in the establishment and maintenance of a condensed chromatin state, also referred to as facultative heterochromatin. BAHD1 is involved in heterochromatin formation and associates with the inactive X chromosome (Xi) in female cells. [PMID:19666599]

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

