

Allgemeine Informationen

Katalog-Nr.:

24000-1-AP

Größe:

150ul, Konzentration: 300 µg/ml von Nanodrop und 267 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG21202

GenBank-Zugangsnummer:

BC028840

GeneID (NCBI):

81573

Vollständiger Name:

ankyrin repeat domain 13C

Berechnete Masse:

541 aa, 61 kDa

Beobachtete Masse:

65-70 kDa

Reinigungsmethode:

Antigen-affinitätsgereinigt

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB : HeLa-Zellen, Maushodengewebe, MDA-MB-231-Zellen

Hintergrundinformationen

ANKRD13C, also termed as ankyrin repeat domain 13C, belongs to the ankyrin repeat domain 13 family. ANKRD13C is localized on ER membranes where it binds to DP to control its biogenesis and trafficking. ANKRD13C modulates the expression and maturation of other G protein coupled receptors as well, indicating that ANKRD13C is a novel regulator of the biogenesis and trafficking of G protein coupled receptors through the biosynthetic pathway.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



HeLa cells were subjected to SDS PAGE followed by western blot with 24000-1-AP (ANKRD13C antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.