

Nur für Forschungszwecke

MRPL47 Polyklonaler Antikörper

Katalog-Nr.: 24728-1-AP



Allgemeine Informationen

Katalog-Nr.:
24728-1-AP

Größe:

150ul, Konzentration: 1000 µg/ml von 57129

Nanodrop und 480 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG11744

GenBank-Zugangsnummer:

BC032522

GeneID (NCBI):

57129

Vollständiger Name:

mitochondrial ribosomal protein L47

Berechnete Masse:

230 aa, 27 kDa

Beobachtete Masse:

29 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB : Maushodengewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



mouse testis tissue were subjected to SDS PAGE followed by western blot with 24728-1-AP (MRPL47 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.