

DTWD2 Polyklonaler Antikörper

Katalog-Nr.: 24873-1-AP

Allgemeine Informationen

Katalog-Nr.:

24873-1-AP

Größe:

150ul, Konzentration: 400 µg/ml von Nanodrop und 320 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG21665

GenBank-Zugangsnummer:

BC114955

GeneID (NCBI):

285605

Vollständiger Name:

DTW domain containing 2

Berechnete Masse:

298 aa, 33 kDa

Beobachtete Masse:

33 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

IP 0.5-4.0 ug für IP und 1:200-1:1000 für WB

Anwendungen

Geprüfte Anwendungen:

IP, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

IP: COLO 320-Zellen,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



IP Result of anti-DTWD2 (IP:24873-1-AP, 3ug;
Detection:24873-1-AP 1:300) with COLO 320 cells
lysate 2000ug.