

Allgemeine Informationen

Katalog-Nr.:

24891-1-AP

Größe:

150ul, Konzentration: 900 µg/ml von Nanodrop und 520 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG19511

GenBank-Zugangsnummer:

BC005926

GeneID (NCBI):

2124

Vollständiger Name:

ecotropic viral integration site 2B

Berechnete Masse:

448 aa, 49 kDa

Beobachtete Masse:

65 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

IHC 1:20-1:200

Anwendungen

Geprüfte Anwendungen:

FC, IHC, WB, ELISA

In Publikationen genannte Anwendungen:

IHC

Getestete Reaktivität:

Human

Zitierte Arten:

Human

Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.

Positivkontrollen:

WB : HepG2-Zellen, A549-Zellen

IHC : humanes Tonsillitisgewebe,

Hintergrundinformationen

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Bingsheng Yang	33574919	Oncol Lett	IHC
Tianyu Xie	36511103	Exp Biol Med (Maywood)	IHC

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




HepG2 cells were subjected to SDS PAGE followed by western blot with 24891-1-AP (EVI2B Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



1X10⁶ Jurkat cells were stained with 0.5ug EVI2B antibody (24891-1-AP, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 24891-1-AP (EVI2B Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 24891-1-AP (EVI2B Antibody) at dilution of 1:50 (under 40x lens).