

Allgemeine Informationen

Katalog-Nr.:
25152-1-AP

Größe:

150ul, Konzentration: 1300 µg/ml von 56242

Nanodrop und 627 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG18570

GenBank-Zugangsnummer:

BC125063

GeneID (NCBI):

56242

Vollständiger Name:

zinc finger protein 253

Berechnete Masse:

499 aa, 58 kDa

Beobachtete Masse:

58-65 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB: HepG2-Zellen, SMMC-7721-Zellen

Hintergrundinformationen

ZNF253, also named as Zinc finger protein 411 or BMZF1, is a 499 amino acid protein, which contains one KRAB domain and eleven C2H2-type zinc fingers. ZNF253 belongs to the krueppel C2H2-type zinc-finger protein family and is expressed in bone marrow and in monocytic and immature erythroid cell lines. ZNF253 may function as a transcription factor.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



HepG2 cells were subjected to SDS PAGE followed by western blot with 25152-1-AP (ZNF253 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.