

Nur für Forschungszwecke

ZNF689 Polyklonaler Antikörper

Katalog-Nr.: 25573-1-AP



Allgemeine Informationen

Katalog-Nr.:

25573-1-AP

Größe:

150ul, Konzentration: 900 µg/ml von Nanodrop und 353 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG22252

GenBank-Zugangsnummer:

BC014000

GeneID (NCBI):

115509

Vollständiger Name:

zinc finger protein 689

Berechnete Masse:

500 aa, 57 kDa

Beobachtete Masse:

50-57 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB: HeLa-Zellen, HepG2-Zellen, Maushodengewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



HeLa cells were subjected to SDS PAGE followed by western blot with 25573-1-AP (ZNF689 Antibody) at dilution of 1:600 incubated at 4 degree celsius over night.