

Nur für Forschungszwecke

SLC12A9 Polyklonaler Antikörper

Katalog-Nr.:26827-1-AP



Allgemeine Informationen

Katalog-Nr.:
26827-1-AP

Größe:
150ul

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG25351

GenBank-Zugangsnummer:
BC000154

GeneID (NCBI):
56996

Vollständiger Name:
solute carrier family 12
(potassium/chloride transporters),
member 9

Berechnete Masse:
96 kDa

Beobachtete Masse:
101 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human, Maus

Positivkontrollen:

WB : HEK-293-Zellen, LNCaP-Zellen,
Mausnierengewebe, MCF-7-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 26827-1-AP (SLC12A9 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.